



SPECIAL APPLICATIONS



**Stainless
steel
version**

s.s AISI 316 for all parts in contact with product

s.s AISI 304 for all the other parts





Dosing Station for liquid/jelly filling



Filling nozzle with "no capsule - no fill" device

Jacketed heated hopper

Dosing adjustment



In-Cap® with 3 pellets dosing stations

1st 2nd 3rd dosing groups for pellets



Powder dosing station with reduced bowl




Upper powder hopper

Lower powder hopper positioned under the dosing disc



Mini-tablets: double filling station for precise counting



a b a Double hopper
Transferring pipes (a)
into filling nozzle (b)

Setting for
precise counting

Control panel