


# CAPPLUS TECHNOLOGIES

## LABORATORY BLOOM VISCOSITY TESTER VISCOSIMETER RVDV-II PRO



CLOSE-UP VIEW OF  
CONTROL PANEL

FEATURES	BENEFITS
CONTINUOUS DISPLAY	SHOWS: VISCOSITY (CP OR MPa-S), TEMPERATURE, SHEAR RATE, SHEAR, STRESS, % TORQUE, SPINDLE
MEETS 4 RECOGNIZED STANDARDS	FOOD ADDITIVE GELATIN GB6783-94, HARD GELATIN CAPSULE FOR MEDICINE GB1371-92, GELATIN FOR MEDICINE QB2354-99, BLOOM VISCOSITY TESTER Q/12XQO193-2004
54 SELECTABLE SPEEDS	PROVIDE SUPERIOR RANGE OF VISCOSITY / SHEAR MEASUREMENTS FOR USE IN STAND-ALONE OPERATION
DV LOADER SOFTWARE (INCLUDED)	ENABLES DOWNLOAD OF CUSTOM SOFTWARE
OPTIONAL WINGATHER™ SOFTWARE	ENABLES AUTOMATIC DATA COLLECTION AND HISTORICAL COMPARISON
BUILT-IN RTD TEMPERATURE PROBE	EASY SAMPLE MONITORING
KEYPAD CONTROL	SIMPLE SELECTION OF TEST PARAMETERS
AUTO-ZERO FUNCTION	ENSURES PRECISION TORQUE MEASUREMENT
MEASURES UP TO 99.9 MPa (MEGAPASCALS)	HANDLES A WIDE VARIETY OF SOLUTION VISCOSITIES
ACCURATE TORQUE MEASUREMENT	+/- 1% OVER SCALE RANGE
REPEATABILITY	+/- 0.2% OVER SCALE RANGE